

Repetitive Pipette HandyStep® S, BRAND



Serial pipetting is quick and easy with the manual HandyStep® S repeating pipette. When used with BRAND PD-Tip syringe tips, the instrument can dispense volumes from 2 µl to 5 ml up to 49 times from a single aspiration.

- Stroke length setting wheel with 9 positions – simple one hand adjustment for right and left-handed person
- 59 different partial volumes can be set with 10 sizes of PD-Tips
- 2-page volume table at the back of the device for quick volume calculation
- Easy connection of tips, the PD-Tips are inserted straight from below
- Safe tip removal at the push of a button; no need to touch the dispenser tips
- Easy-to-read display, no battery required
- Innovative plastics for increased resistance to chemicals
- Maintenance-free, factory-adjusted
- Best results when used with BRAND PD-Tip syringe tips. Combitips® plus, Combitips®, Repet-Tips, Encode™ tips and other compatible pipette tips can be used as well.

Ideal for routine or research applications in diagnostics, molecular biology, environmental analysis, and many more. HandyStep® S and PD-tips work on the direct displacement principle. This offers the highest-precision dispensing of liquid media with high viscosity, high density, or high vapor pressure.

Items supplied:

HandyStep® S, DE-M marked, with serial number, performance certificate, operating manual, shelf/rack mount, 3 PD-Tips: 0.1 ml, 1 ml and 10 ml.

Article	Order No.
① HandyStep®, incl. shelf/rack mount	10 49 75110
Shelf/Rack mount for HandyStep® S	10 49 75130

